

Vital signal Detection with earphone(2)

1. To detect the vital signal from anyone's ear

- (a) Basically almost canal type inside earphones are possible to use
- (b) For the open-ear type earphone, an adapter into canal type
- (c) Almost headphones are possible to use
Over the hair, some of them are possible to use
- (d) Different sizes from many companies are helpful





But for the people who has still difficulty


2. Salustek has decided to adopt the 3D scanner from United Sciences LLC




During scanning



Entire scan image



Pick the necessary part up and modify if necessary and make stl file



With 3D printing from stl file



Equip a speaker unit



Salustek will start its service for the person who has difficulty to detect vital signals with his or her own earphone.

We will announce the detail on Salustek home page.